

Ambu® UltraSeal

DATASHEET

Ambu UltraSeal is disposable face mask designed specifically for anesthetic departments.

The anatomically shaped cuff secures a tight seal and enables optimal ventilation. The face mask has good edge grip which makes it easy to use while providing a tight seal. The mask shape ensures that the inner shape of the dome is spacious enough for the mouth and nose also when pressure is applied to the mask.

The UltraSeal face mask is not made with natural rubber latex and is MR safe.

Complete line with 7 different sizes from size 0 Pediatric to size 6 Large adult. All sizes come with and without check valve and some sizes are available with hook ring.



KEY BENEFITS

Very soft anatomically shaped cuff enabling a tight seal with minimum applied pressure

Flexible, self-expanding dome

Shoulder grip fitting different hand sizes

Crystal clear dome for easy observation of the patient's condition

The single patient use concept offers minimized risk of cross infection

All sizes are delivered individually packed in a transparent, easy to open bag

Color coding on mask label for fast and easy identification of size

All sizes are delivered individually packed, non-sterile, in packs of 20 units

Phtalate-free material

Available with or without the hook ring

MR safe

Ambu

ORDERING SPECIFICATIONS WITHOUT HOOK RING

Item no.	Packaging Color Code	Type
305102000	●	Ambu UltraSeal # 0 Pediatric
305104000	●	Ambu UltraSeal # 1 Pediatric
305106000	●	Ambu UltraSeal # 2 Pediatric
305108000H	●	Ambu UltraSeal # 3 Pediatric
305010000H	●	Ambu UltraSeal # 4 Small adult/Child
305012000H	●	Ambu UltraSeal # 5 Medium adult
305014000H	●	Ambu UltraSeal # 6 Large adult
305101000	●	Ambu UltraSeal # 0 Pediatric, with check valve
305103000	●	Ambu UltraSeal # 1 Pediatric, with check valve
305105000	●	Ambu UltraSeal # 2 Pediatric, with check valve
305107000H	●	Ambu UltraSeal # 3 Pediatric, with check valve
305011000H	●	Ambu UltraSeal # 4 Small adult/Child, with check valve
305013000H	●	Ambu UltraSeal # 5 Medium adult, with check valve
305015000H	●	Ambu UltraSeal # 6 Large adult, with check valve

WITH HOOK RING

Item no.	Packaging Color Code	Type
305108000	●	Ambu UltraSeal # 3 Pediatric
305010000	●	Ambu UltraSeal # 4 Small adult/Child
305012000	●	Ambu UltraSeal # 5 Medium adult
305014000	●	Ambu UltraSeal # 6 Large adult
305107000	●	Ambu UltraSeal # 3 Pediatric, with check valve
305011000	●	Ambu UltraSeal # 4 Small adult/Child, with check valve
305013000	●	Ambu UltraSeal # 5 Medium adult, with check valve
305015000	●	Ambu UltraSeal # 6 Large adult, with check valve

SIZE

Type	Weight without check valve	Weight with check valve
# 0 Pediatric	9 g	10 g
# 1 Pediatric	12 g	13 g
# 2 Pediatric	19 g	20 g
# 3 Pediatric	25 g	26 g
# 4 Small adult/Child	44 g	45 g
# 5 Medium adult	53 g	54 g
# 6 Large adult	72 g	73 g

MATERIALS

Product	
Dome	Polyvinylchloride
Cuff	Polyvinylchloride
Check valve	Polycarbonate
Packaging	
Pouch	PP (Polypropylene)
Box	Cardboard
Environment	
Phthalate-free	Yes
Not made with natural rubber latex	Yes



Ambu UltraSeal without check valve, with hook ring.



Ambu UltraSeal with check valve, without hook ring.

Ambu

Ambu A/S

Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com