

Ambu® WhiteSensor 4500M-H

ECG Electrode - Single use



Key Benefits

- Highly conductive wet gel
- Good adhesion
- Flexible foam backing

Ambu® WhiteSensor 4500M-H

The Ambu WhiteSensor 4500M-H features a wet gel with good adhesion to ensure a good signal quality during short-term ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode ensures ease of use and comfort during application.



Wet gel



Foam backing

WhiteSensor

Ambu^{*}
Ideas that work for life

Specifications

Dimensions	
Electrode size (length x width in mm)	50 x 48
Adhesive area (in mm ²)	1585
Height excluding connector (in mm)	1
Sensor	
Sensor material	Silver, silver chloride (Ag/AgCl)
Gel system	Wet gel
Gel area (in mm ²)	201
Sensor area (in mm ²)	79
Electrical data (ANSI/AAMI)	
AC impedance - typical	65 Ω
DC offset potential - typical	1 mV
Defibrillation overload recovery - typical	12 mV
Rate of change for polarisation potential - typical	0.3 mV/s
Combined offset instability and internal noise	9 μV
Bias current tolerance	5 mV

Environment	
PVC-free electrode	Yes
Electrode not made with natural rubber latex	Yes
PVC-free packaging	Yes
X-Ray & MRI	
Radiolucent	No
MR Conditional	No
Shelf life	
Opened pouch	30 days
Unopened pouch	24 months*

* From date of production

Materials

Electrode	
Bio-compatible	Yes
Sponge	Polyester Polyurethane Reticulated
Backing material	Polyethylene foam (PE)
Backing material adhesive	Polyacrylate
Supporting label	Polyethylene
Connector	Stainless Steel
Release liner	Polyester (PET)

Packaging	
Pouch (outer/middle/inner layers)	Polyester/Aluminum/Polythene
Boxes	Corrugated cardboard
Options	
Connectors	S
Precaution	
Single use only	
Accessories	
Skin Preps	Ref: 2121M

Available configurations

Reference	Units/Card	Units/pouch	Units/inner box	Units/outer box
4500M-H	3	60	300	3000
4500M-H/10	5	10	300	3000

Distributor in EU
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 ambu.com

Manufacturer
MEDICO ELECTRODES INTERNATIONAL LTD.
 [AN ISO13485 COMPANY]
 Plot 142A/11, 12, 27, 28 & 29,
 Noida Special Economic Zone,
 Noida, Uttar Pradesh - 201305, India
 www.medicoelectrodes.com

CE US: Rx only