

DATASHEET

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope™ Duodeno, you have a new scope with consistent performance every time, and there is no risk of patient cross-contamination. What's more, there's no need for costly reprocessing or repair.



Distal tip showing elevator in use with a sphincterotome



Elevator control lever to raise (up) and lower (down) the elevator

Ambu® aScope™ Duodeno Sterile. Single-use. Problem solved.



KEY BENEFITS

Sterile

No risk of patient cross-contamination

Cost transparency

Eliminates the risk and cost of reprocessing

Familiar design and function to that of conventional duodenoscopes

Ambu

SPECIFICATIONS

Optical System

| | |
|-------------------|-----------------------|
| Direction of View | Retro Side Viewing 6° |
| Field of View | 130° |

Insertion Section

| | |
|-------------------------------|-----------|
| Distal End Outer Diameter | φ 13.7 mm |
| Insertion Tube Outer Diameter | φ 11.3 mm |
| Working Length | 1240 mm |

Instrument Working Channel

| | |
|----------------|--------------------|
| Inner Diameter | φ 4.2 mm (12.6 FR) |
|----------------|--------------------|

Bending Section

| | |
|------------------|---|
| Angulation Range | Up 120° / Down 90° Left 90° / Right 110° |
|------------------|---|

Illumination

| | |
|---------------------|-----------|
| Illumination source | Dual LEDs |
|---------------------|-----------|

Dimensions


| | |
|--------|-------|
| Weight | 700 g |
|--------|-------|

ORDERING SPECIFICATIONS*

| Reference number | Description | Quantity |
|------------------|----------------|------------|
| 485001000 | aBox Duodeno | 1 |
| 482001000/4 | aScope Duodeno | Box of 1** |
| 482001000 | aScope Duodeno | Box of 6 |

*Some configurations may not be available worldwide

**Each inner box contains one scope and IFU



Ambu

Ambu A/S
 Baltorpbakken 13
 2750 Ballerup
 Denmark
 T +45 72 25 20 00
 ambu.com


 2797

US: Rx Only