

# Ambu<sup>®</sup> Disposable Pressure Manometer



## Key Benefits

- The slim design makes it easy to monitor pressure with minimum obstruction of the view of the patients face
- The single patient use concept offers minimized risk of cross infection
- All materials are Latex and PVC free
- The manometer has a colour coded scale to give enhanced visual support during ventilation. The colour coding scale is only for visual support. The correct ventilation pressure must be determined by the medical professional

The Ambu Disposable Pressure Manometer for single patient use is intended to be used for monitoring the patient's airway pressure.

The Ambu Disposable Pressure Manometer is suitable for use with Ambu resuscitators or other resuscitators, hyperinflation bags, CPAP masks or circuits.

## Specifications

Size	
Length	55 mm (2.17")
Diameter	22 mm (0.87")
Weight	6.4 g
Connector size	ID 3.7 mm (.15")
Materials	
Connector	SEBS
Spring	CuSn
Housing	SB
Piston	PP
Membrane	Silicone
Operating temperature range	-18°C a 50°C
Storage temperature range	-40°C a 60°C
Measuring	
Measuring Limits	0-60 cm H <sub>2</sub> O (hPa)
Accuracy	+/- 2 cm H <sub>2</sub> O (hPa) a 7.5, 10, 15, 20, 30 H <sub>2</sub> O (hPa) +/- 3 cm H <sub>2</sub> O (hPa) a 40 H <sub>2</sub> O (hPa) +/- 5 cm H <sub>2</sub> O (hPa) a 60 H <sub>2</sub> O (hPa)



Ambu Disposable Pressure Manometer connected to Ambu SPUR II Resuscitator

## Ordering specifications

Item no.	
322 003 000	Ambu Disposable Pressure Manometer EU 20 pcs
322 005 000	Ambu Disposable Pressure Manometer EU 240 pcs.
322 004 000	Ambu Disposable Pressure Manometer US 20 pcs.
322 006 000	Ambu Disposable Pressure Manometer US 240 pcs.

### Ambu Australia Pty Ltd

#### Customer Services

Tel.: 1300 233 118  
Fax. 1300 156 526  
aussales@ambu.com

### Ambu Asia (Sales office)

#### Customer Services

Tel.: +604 252 9100  
Fax. +604 644 8932  
csasia@ambu.com