

Ambu® Man W + I



Key Benefits

- Wireless connections
- Ambu's unique hygienic system
- The rigidity of the chest is adjustable
- Airways open only when head is correctly extended
- Jaw thrust possible
- PC compatible
- Training mat integrated in the carrying bag

Ambu® Man W + I

The AmbuMan is a highly developed instruction- and training manikin for the simulation of realistic and correct anatomical and physiological conditions during cardiopulmonary resuscitation by means of rescue breathing and chest compression.

Realistic scenarios

The AmbuMan W (Wireless) does not need cables and wires. Which means it can be placed in realistic surroundings away from computers and power supplies. This gives a great amount of freedom during training and because there are no wires or cables lying around it gives an even more realistic scenario.

Low risk of cross-contamination

The AmbuMan incorporates an unique, patented hygienic system under which all trainees get their own face piece and head bag, thereby minimizing any risk of cross-infection between the trainees.

The hygienic system also features realistic expiration from the AmbuMan through the mouth and nose.

Mechanical Instrument

Version W:

The AmbuMan W models are equipped with a mechanical monitoring instrument as well as digital sensors showing effectiveness of the resuscitation. It also shows results achieved for insufflations volume and depth of external chest compression, as well as it indicates any stomach inflation and incorrect hand positioning.

AmbuMan W models can be connected to a computer through USB or WLAN using the Ambu CPR-Software for a more detailed analysis of the performed resuscitation.

Version I:

The AmbuMan I (Instrument) models are equipped with a mechanical monitoring instrument showing effectiveness of the resuscitation, including results achieved for ventilation volume and depth of external chest compression. Moreover it indicates any stomach inflation and incorrect hand positioning.

Specifications

Item no.	Description
Weight: torso with carrying case	12 kg
Full body with carrying case	17 kg
Total length: Torso	80 cm
Full body	170 cm

Materials

Product	Material
Basic unit	Polyethylene
Instrument part	ABS plastic
Skull	PVC, hard
Chest, face	PVC, soft

Ordering specifications

Item no.	Description
A234 204 000	AmbuMan Model I, Torso
A234 307 000	AmbuMan Model W, Torso with CPR-Software
234 000 706	Upgrade Kit - Fill Body
234 000 704	Upgrade Kit -Arms
Accessories	
234 000 703	5 Face pieces
234 000 702	100 Head bags
234 000 790	Jogging suit jacket
234 000 792	Jogging suit trousers
234 000 786	Legs with trousers
234 000 794	Carrying bag



AmbuMan computer interface



Arms and legs for AmbuMan



Bag for AmbuMan

Ambu Australia Pty Ltd

Customer Services

Tel.: 1300 233 118

Fax: 1300 156 526

aussales@ambu.com

Ambu Asia (Sales office)

Customer Services

Tel.: +604 252 9100

Fax: +604 644 8932

csasia@ambu.com



Contains FCC ID U30-G2M5477